

BJOY RING

Access that travels with you

# Joaquín, our first user

*“Before, it was enough for me to access a computer. Now I want to access different devices and I need to do it efficiently and comfortably.”*



# The idea that inspired us

---

If a joystick is already a  
convenient means of  
access...

Why change?




# Installation

- ❑ **Self-installable**
- ❑ **Compatible with the different joysticks**
- ❑ **Does not interfere** with the wheelchair's electronics (based on accelerometers)




# Personalisation



- BJOY Ring can be configured on start-up and all the parameters are **saved to the device**.
- The device is **calibrated** and adjusts to the user's position and movements.

# Personalisation



- **Clicking** can be done with external **switches** or by **dwell click** (without requiring software).
- Provides **ways** for double-clicking, right-button clicking, drag&drop, scroll, etc.

# Two versions available

## BJOY Ring (wired)

- Ideal for controlling devices installed in the wheelchair itself.



## BJOY Ring Wireless

- Ideal for controlling elements that are not in the wheelchair or many devices at the same time.



# BJOY Ring Wireless. Multiple receivers

---



Just getting close is enough



# BJOY Ring Wireless. Multiple receivers

## □ Scenario 2: Nearby devices



Toggle between devices with a single **gesture** (making a circle with the wheelchair's joystick) or clicking on the external switch.

# BJOY Ring Wireless. Multiple receivers

- ❑ **Scenarios 1 and 2 can be combined for up to a total of four devices.**



# Technical Specifications

## General

- Plug and play USB 2.0
- Compatible with Windows, Linux, Chrome OS, Mac OS and Android
- Automatic angle recognition
- Initial point between 0 and 45°
- Two switch inputs
- Acoustic feedback

## Wired version


- Extension USB cable included

## Wireless version

- Up to 4 receivers per emitter
- 2 metres of coverage
- Rechargeable battery (micro USB)
- Use it also as the wired version



# What our users are saying...



*"It comes as natural to me as manoeuvring my wheelchair. It helps me study"*

*"We were surprised by how easy it was to install"*

*"Alternating between my mobile and computer control is really convenient"*

*"I have gained autonomy. No one has to help me"*

*"I can now use my social networks from my wheelchair, so I'm connected to my people and the things I'm interested in"*

“GIVE ME ACCESS TO THE  
INTERNET AND I’LL CHANGE THE  
WORLD”

Thanks for your attention. Find us at **[www.bjliveat.com](http://www.bjliveat.com)**